RESMED



One fit. One solution. One for all.

The perfect balance of sophisticated design and easy functionality, Mirage™ FX is the solution you've been waiting for: performance, ease-of-use and quality all rolled into one.

Benefits for you ...

Sizing is a one-step process, with Mirage FX covering a wide range of facial characteristics in one mask ... so whoever your patients are, Mirage FX is the one for all

Fitting is just as effortless. Mirage FX is intuitive to fit and adjust, so fitting takes virtually no time at all. With only four parts, including headgear, you'll also save time on managing inventory.

Mirage FX—puts your business first.

... and your patients

The Mirage FX Spring Air™ cushion provides even pressure distribution and comfort for your patients, incorporating extra softness at the nasal bridge to reduce discomfort commonly felt in this area.

Your patients will also find Mirage FX easy to live with: it's easy to clean, wear and adjust ... and being one of the lightest nasal masks on the market, it won't weigh them down.

Mirage FX—make it your number one mask.



RESMED

Mirage FX is everything you've been waiting for all rolled into one.



effortless fit

Intuitive adjustments save time and mean just about anyone can achieve a perfect, first-time fit.

simple solution

Mirage FX is the solution you've been waiting for: performance, ease-of-use and quality all rolled into one.

mask for all

With Mirage FX, you'll fit more people the first time round with one mask, saving time, money and additional patient visits.



New, form-fitting headgear springs to life, ready to fit and wear



Flexible, winged forehead support allows a dynamic range of fit and adjustment



SoftEdge™ headgear design is kind to the face, reducing marks and discomfort

Lightweight, durable frame

strong enough to stabilise the mask on the face, soft enough to accommodate user movement during sleep



provides even pressure distribution, with extra softness at the nasal bridge, reducing common discomfort in this sensitive area

Diffuse vent design

disperses air gently and quietly, away from bed partner

Headgear loops

simplify attachment and detachment of headgear; retain optimal headgear settings

Squeeze-tab elbow

makes it easy to attach and detach the elbow from the frame

Complete System

	Europe 1	Europe 2	Europe 3
Mirage FX	62100	62101	62102

Europe 1: English (UK), French, German and Italian (Switzerland) Europe 2: English, Danish, Finnish, Norwegian and Swedish

Europe 3: English, Czech, Dutch (Belgium), French (Belgium), German (Belgium), Greek, Hungarian, Italian, Polish, Portuguese, Russian, Slovak, Spanish and Turkish Mirage[™] FX NASAL MASK

www.resmedmiragefx.com

Manufacturer: ResMed Ltd 1 Elizabeth Macarthur Drive Bella Vista NSW 2153 Australia. Distributed by: ResMed Corp 9001 Spectrum Center Boulevard San Diego CA 92123 USA. ResMed (UK) Ltd 96 Milton Park Abingdon Oxfordshire OX14 4RY UK. See www.resmed.com for other ResMed locations world.wide. Mirage and Spring Air are trademarks of ResMed Ltd and are registered in U.S. Patent and Trademark Office. © 2010 ResMed Ltd. Specifications may change without notice. 1014062/3 2010-11